

Pulseq Tutorial Notes


(01_from_FID_to_PRESS)

Qingping Chen

*Division of Medical Physics, Dept. Of Radiology,
University Medical Center Freiburg, Germany*

June 20 2023

Hard RF pulse – mr.makeBlockPulse()




$BW = 1 / (4 * Duration)$

$delay = \max(user_desired_delay, deadTime)$

$event\ duration = RF_delay + shape_dur + ringdownTime$

ADC – mr.makeAdc()




$\text{delay} = \max(\text{user_defined_delay}, \text{deadTime})$

$\text{event duration} = \text{delay} + \text{num_sample} * \text{ADCdwellTime} + \text{deadTime}$

The number of ADC samples (num_sample) on Siemens need to be divisible by 4.

Note that both Pulseseq and Siemens define the ADC samples to happen in the centre of the dwell periods.


Trapezoid gradient – mr.makeTrapezoid()



`grad.area=amplitude*(flatTime+riseTime/2+fallTime/2)`


`grad.flatArea=amplitude*flatTime`

Sinc RF pulse – mr.makeSincPulse()



```
BW=timeBwProduct/duration
alpha=apodization
N=round(duration/dwell)
t=(1:N - 0.5)*dwell
tt=t-duration*centerpos
window=(1-alpha+alpha*cos(2*pi*tt/duration))
signal=window*sinc(BW*tt)
flip=sum(signal)*dwell*2*pi
signal=signal*flipAngle/flip
amplitude=BW/sliceThickness
area=amplitude*duration
gz_flatArea=area
gz_flatTime=duration
rf_delay=gz_delay+gz_riseTime
gxr_area=area*(1-centerpos)-0.5*gz_area-area
```

Extended trapezoid gradient – mr.makeExtendedTrapezoid()



Example: slice-selective gradient with crushers, here composed from 3 gradients:

gradient_1

Amplitudes: [a, b, c, d]

Times: [t_a, t_b, t_c, t_d]

gradient_2

Amplitudes: [d, e]

Times: [t_d, t_e]

gradient_3


Amplitudes: [e, f, g, h]

Times: [t_e, t_f, t_g, t_h]


Slice-selective excitation with crushers

`makeExtendedTrapezoid()` creates an “extended trapezoid” gradient event by specifying a set of points (amplitudes) at specified time points (times) at a given channel with given system limits. The amplitude at the beginning and the end do not have to be 0.


s01_FID (single flip angle)




s02_FID_multipleFAs




s03_SE



s04_SE_withBackGrad



s05_withSpoilers



s06_PRESS

